

Antibodies Product Datasheet

Product Name:	PAR1 (Cleaved-Ser42) Antibody
Catalog Number:	L0374sec
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human PAR1 (Cleaved-Ser42) Polyclonal
Concentration:	1 mg/ml
Gene Name:	F2R CF2R PAR1 TR
Protein Name:	Proteinase-activated receptor 1 (PAR-1) (Coagulation factor II receptor) (Thrombin receptor)
Other Name:	Proteinase-activated receptor 1 (PAR-1) (Coagulation factor II receptor) (Thrombin receptor)
Human Swiss Prot No.:	P25116
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PAR1 (Cleaved-Ser42).
Storage and Stability:	-20°C/1 year