

Antibodies Product Datasheet

Product Name:	Collagen V α 1 (Cleaved-Ala1605) Antibody
Catalog Number:	L0386
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Collagen V α 1 (Cleaved-Ala1605)
Concentration:	1 mg/ml
Gene Name:	COL5A1
Protein Name:	Collagen V α 1 (Cleaved-Ala1605)
Other Name:	Collagen alpha-1(V) chain
Human Swiss Prot No.:	P20908
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Collagen V α 1 (Cleaved-Ala1605, protein was cleaved amino acid sequence between 1605-1606)
Storage and Stability:	-20°C/1 year