

Antibodies Product Datasheet

Product Name:	ITGA5 (heavy chain, Cleaved-Arg894) Antibody
Catalog Number:	L0392
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ITGA5 (heavy chain, Cleaved-Arg894)
Concentration:	1 mg/ml
Gene Name:	ITGA5 FNRA
Protein Name:	ITGA5 (heavy chain, Cleaved-Arg894)
Other Name:	Integrin alpha-5 (CD49 antigen-like family member E) (Fibronectin receptor subunit alpha) (Integrin alpha-F) (VLA-5) (CD antigen CD49e) [Cleaved into: Integrin alpha-5 heavy chain; Integrin alpha-5 light chain]
Human Swiss Prot No.:	P08648
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse ITGA5 (heavy chain, Cleaved-Arg894, protein was cleaved amino acid sequence between 894-895)
Storage and Stability:	-20°C/1 year