

Antibodies Product Datasheet

Product Name:	ITGAE (light chain, Cleaved-Arg177) Antibody
Catalog Number:	L0395
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ITGAE (light chain, Cleaved-Arg177)
Concentration:	1 mg/ml
Gene Name:	ITGAE
Protein Name:	ITGAE (light chain, Cleaved-Arg177)
Other Name:	Integrin alpha-E (HML-1 antigen) (Integrin alpha-IEL) (Mucosal lymphocyte 1 antigen) (CD antigen CD103) [Cleaved into: Integrin alpha-E light chain; Integrin alpha-E heavy chain]
Human Swiss Prot No.:	P38570
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human ITGAE (light chain, Cleaved-Arg177, protein was cleaved amino acid sequence between 177-178)
Storage and Stability:	-20°C/1 year