

Antibodies Product Datasheet

Product Name:	ITGAV (light chain, Cleaved-Asp891) Antibody
Catalog Number:	L0396
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ITGAV (light chain, Cleaved-Asp891)
Concentration:	1 mg/ml
Gene Name:	ITGAV MSK8 VNRA
Protein Name:	ITGAV (light chain, Cleaved-Asp891)
Other Name:	Integrin alpha-V (Vitronectin receptor subunit alpha) (CD antigen CD51) [Cleaved into: Integrin alpha-V heavy chain; Integrin alpha-V light chain]
Human Swiss Prot No.:	P06756
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse ITGAV (light chain, Cleaved-Asp891, protein was cleaved amino acid sequence between 890-891)
Storage and Stability:	-20°C/1 year