

Antibodies Product Datasheet

Product Name:	CFAD (Cleaved-Ile26) Antibody
Catalog Number:	L0418
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CFAD (Cleaved-Ile26)
Concentration:	1 mg/ml
Gene Name:	CFD DF PFD
Protein Name:	CFAD (Cleaved-Ile26)
Other Name:	Complement factor D (EC 3.4.21.46) (Adipsin) (C3 convertase activator) (Properdin factor D)
Human Swiss Prot No.:	P00746
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human CFAD (Cleaved-Ile26, protein was cleaved amino acid sequence between 25-26)
Storage and Stability:	-20°C/1 year