

Antibodies Product Datasheet

Product Name:	Eosinophil peroxidase (heavy chain, Cleaved-Val251) Antibody
Catalog Number:	L0421
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Eosinophil peroxidase (heavy chain, Cleaved-Val251)
Concentration:	1 mg/ml
Gene Name:	EPX EPER EPO EPP
Protein Name:	Eosinophil peroxidase (heavy chain, Cleaved-Val251)
Other Name:	Eosinophil peroxidase (EPO) (EC 1.11.1.7) [Cleaved into: Eosinophil peroxidase light chain; Eosinophil peroxidase heavy chain]
Human Swiss Prot No.:	P11678
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat Eosinophil peroxidase (heavy chain, Cleaved-Val251, protein was cleaved amino acid sequence between 250-251)
Storage and Stability:	-20°C/1 year