

Antibodies Product Datasheet

Product Name:	FA9 (heavy chain, Cleaved-Val227) Antibody
Catalog Number:	L0427
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FA9 (heavy chain, Cleaved-Val227)
Concentration:	1 mg/ml
Gene Name:	F9
Protein Name:	FA9 (heavy chain, Cleaved-Val227)
Other Name:	Coagulation factor IX (EC 3.4.21.22) (Christmas factor) (Plasma thromboplastin component) (PTC) [Cleaved into: Coagulation factor IXa light chain; Coagulation factor IXa heavy chain]
Human Swiss Prot No.:	P00740
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat FA9 (heavy chain, Cleaved-Val227, protein was cleaved amino acid sequence between 226-227)
Storage and Stability:	-20°C/1 year