

Antibodies Product Datasheet

Product Name:	ITGA8 (heavy chain, Cleaved-Arg906) Antibody
Catalog Number:	L0431
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ITGA8 (heavy chain, Cleaved-Arg906)
Concentration:	1 mg/ml
Gene Name:	ITGA8
Protein Name:	ITGA8 (heavy chain, Cleaved-Arg906)
Other Name:	Integrin alpha-8 [Cleaved into: Integrin alpha-8 heavy chain; Integrin alpha-8 light chain]
Human Swiss Prot No.:	P53708
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat ITGA8 (heavy chain, Cleaved-Arg906, protein was cleaved amino acid sequence between 906-907)
Storage and Stability:	-20°C/1 year