

Antibodies Product Datasheet

Product Name:	Kallikrein-4 (Cleaved-Ile31) Antibody
Catalog Number:	L0434
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Kallikrein-4 (Cleaved-Ile31)
Concentration:	1 mg/ml
Gene Name:	KLK4 EMSP1 PRSS17 PSTS
Protein Name:	Kallikrein-4 (Cleaved-Ile31)
Other Name:	Kallikrein-4 (EC 3.4.21.-) (Enamel matrix serine proteinase 1) (Kallikrein-like protein 1) (KLK-L1) (Prostase) (Serine protease 17)
Human Swiss Prot No.:	Q9Y5K2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Kallikrein-4 (Cleaved-Ile31, protein was cleaved amino acid sequence between 30-31)
Storage and Stability:	-20°C/1 year