

Antibodies Product Datasheet

Product Name:	Kallikrein-6 (Cleaved-Leu22) Antibody
Catalog Number:	L0435
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Kallikrein-6 (Cleaved-Leu22)
Concentration:	1 mg/ml
Gene Name:	KLK6 PRSS18 PRSS9
Protein Name:	Kallikrein-6 (Cleaved-Leu22)
Other Name:	Kallikrein-6 (EC 3.4.21.-) (Neurosin) (Protease M) (SP59) (Serine protease 18) (Serine protease 9) (Zyme)
Human Swiss Prot No.:	Q92876
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Kallikrein-6 (Cleaved-Leu22, protein was cleaved amino acid sequence between 21-22)
Storage and Stability:	-20°C/1 year