

Antibodies Product Datasheet

Product Name:	Kallikrein-7 (Cleaved-Ile30) Antibody
Catalog Number:	L0436
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Kallikrein-7 (Cleaved-Ile30)
Concentration:	1 mg/ml
Gene Name:	KLK7 PRSS6 SCCE
Protein Name:	Kallikrein-7 (Cleaved-Ile30)
Other Name:	Kallikrein-7 (hK7) (EC 3.4.21.117) (Serine protease 6) (Stratum corneum chymotryptic enzyme) (hSCCE)
Human Swiss Prot No.:	P49862
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat Kallikrein-7 (Cleaved-Ile30, protein was cleaved amino acid sequence between 29-30)
Storage and Stability:	-20°C/1 year