

Antibodies Product Datasheet

Product Name:	KLK1/2/3 (Cleaved-Ile25) Antibody
Catalog Number:	L0437
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KLK1/2/3 (Cleaved-Ile25)
Concentration:	1 mg/ml
Gene Name:	KLK1/2/3
Protein Name:	KLK1/2/3 (Cleaved-Ile25)
Other Name:	Kallikrein-2 (EC 3.4.21.35) (Glandular kallikrein-1) (hGK-1) (Tissue kallikrein-2)
Human Swiss Prot No.:	P20151/P06870/P07288
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human KLK1/2/3 (Cleaved-Ile25, protein was cleaved amino acid sequence between 24-25)
Storage and Stability:	-20°C/1 year