

Antibodies Product Datasheet

Product Name:	FA10 (light chain, Cleaved-Arg179) Antibody
Catalog Number:	L0446
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FA10 (light chain, Cleaved-Arg179)
Concentration:	1 mg/ml
Gene Name:	F10
Protein Name:	FA10 (light chain, Cleaved-Arg179)
Other Name:	Coagulation factor X (EC 3.4.21.6) (Stuart factor) (Stuart-Prower factor) [Cleaved into: Factor X light chain; Factor X heavy chain; Activated factor Xa heavy chain]
Human Swiss Prot No.:	P00742
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human FA10 (light chain, Cleaved-Arg179, protein was cleaved amino acid sequence between)
Storage and Stability:	-20°C/1 year