

Antibodies Product Datasheet

Product Name:	FA7 (light chain, Cleaved-Ala61) Antibody
Catalog Number:	L0447
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FA7 (light chain, Cleaved-Ala61)
Concentration:	1 mg/ml
Gene Name:	F7
Protein Name:	FA7 (light chain, Cleaved-Ala61)
Other Name:	Coagulation factor VII (EC 3.4.21.21) (Proconvertin) (Serum prothrombin conversion accelerator) (SPCA) (Eptacog alfa) [Cleaved into: Factor VII light chain; Factor VII heavy chain]
Human Swiss Prot No.:	P08709
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human FA7 (light chain, Cleaved-Ala61, protein was cleaved amino acid sequence between 60-61)
Storage and Stability:	-20°C/1 year