

Antibodies Product Datasheet

Product Name:	FA8A (heavy chain 200k/92k, Cleaved-Ala20) Antibody
Catalog Number:	L0448
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FA8A (heavy chain 200k/92k, Cleaved-Ala20)
Concentration:	1 mg/ml
Gene Name:	F8 F8C
Protein Name:	FA8A (heavy chain 200k/92k, Cleaved-Ala20)
Other Name:	Coagulation factor VIII (Antihemophilic factor) (AHF) (Procoagulant component) [Cleaved into: Factor VIIIa heavy chain, 200 kDa isoform; Factor VIIIa heavy chain, 92 kDa isoform; Factor VIII B chain; Factor VIIIa light chain]
Human Swiss Prot No.:	P00451
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse FA8A (heavy chain 200k or 92k, Cleaved-Ala20, protein was cleaved amino acid sequence between 19-20)
Storage and Stability:	-20°C/1 year