

Antibodies Product Datasheet

Product Name:	CATL1 (heavy chain, Cleaved-Ala114) Antibody
Catalog Number:	L0459
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CATL1 (heavy chain, Cleaved-Ala114)
Concentration:	1 mg/ml
Gene Name:	CTSL1 CTSL
Protein Name:	CATL1 (heavy chain, Cleaved-Ala114)
Other Name:	Cathepsin L1 (EC 3.4.22.15) (Major excreted protein) (MEP) [Cleaved into: Cathepsin L1 heavy chain; Cathepsin L1 light chain]
Human Swiss Prot No.:	P07711
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human CATL1 (heavy chain, Cleaved-Ala114, protein was cleaved amino acid sequence between 113-114)
Storage and Stability:	-20°C/1 year