

Antibodies Product Datasheet

Product Name:	ENTK (light chain, Cleaved-Ile785) Antibody
Catalog Number:	L0462
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ENTK (light chain, Cleaved-Ile785)
Concentration:	1 mg/ml
Gene Name:	TMPRSS15 ENTK PRSS7
Protein Name:	ENTK (light chain, Cleaved-Ile785)
Other Name:	Enteropeptidase (EC 3.4.21.9) (Enterokinase) (Serine protease 7) (Transmembrane protease serine 15) [Cleaved into: Enteropeptidase non-catalytic heavy chain; Enteropeptidase catalytic light chain]
Human Swiss Prot No.:	P98073
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human ENTK (light chain, Cleaved-Ile785, protein was cleaved amino acid sequence between 785-786)
Storage and Stability:	-20°C/1 year