

Antibodies Product Datasheet

Product Name:	NF- κ B p105 (Cleaved-Thr434) Antibody
Catalog Number:	L0473
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human NF- κ B p105 (Cleaved-Thr434)
Concentration:	1 mg/ml
Gene Name:	NFKB1
Protein Name:	NF- κ B p105 (Cleaved-Thr434)
Other Name:	Nuclear factor NF-kappa-B p105 subunit (DNA-binding factor KBF1) (EBP-1) (Nuclear factor of kappa light polypeptide gene enhancer in B-cells 1) [Cleaved into: Nuclear factor NF-kappa-B p50 subunit]
Human Swiss Prot No.:	P19838
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse NF- κ B p105 (Cleaved-Thr434, protein was cleaved amino acid sequence between 433-434)
Storage and Stability:	-20°C/1 year