

Antibodies Product Datasheet

Product Name:	Caspase 4 (p10, Cleaved-Ala290) Antibody
Catalog Number:	L0480
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Caspase 4 (p10, Cleaved-Ala290)
Concentration:	1 mg/ml
Gene Name:	CASP4 ICH2
Protein Name:	Caspase4
Other Name:	Caspase-4 (CASP-4) (EC 3.4.22.57) (ICE(rel)-II) (Protease ICH-2) (Protease TX) [Cleaved into: Caspase-4 subunit 1; Caspase-4 subunit 2]
Human Swiss Prot No.:	P49662
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Caspase 4 (p10, Cleaved-Ala290, protein was cleaved amino acid sequence between 289-290)
Storage and Stability:	-20°C/1 year