

Antibodies Product Datasheet

Product Name:	ApoA-I Monoclonal Antibody
Catalog Number:	N0030
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human ApoA-I expressed in E. Coli.
Concentration:	
Gene Name:	APOA1
Protein Name:	Apolipoprotein A-I
Other Name:	APOA1; Apolipoprotein A-I; Apo-AI; ApoA-I; Apolipoprotein A1
Human Swiss Prot No.:	P02647
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	ApoA-I Monoclonal Antibody detects endogenous levels of ApoA-I protein.
Storage and Stability:	-20°C/1 year