

Antibodies Product Datasheet

Product Name:	APP Monoclonal Antibody
Catalog Number:	N0041
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of APP expressed in E. Coli.
Concentration:	
Gene Name:	APP
Protein Name:	Amyloid beta A4 protein, Amyloid- β , A β
Other Name:	APP; A4; AD1; Amyloid beta A4 protein; ABPP; APPI; APP; Alzheimer disease amyloid protein; Cerebral vascular amyloid peptide; CVAP; PreA4; Protease nexin-II; PN-II
Human Swiss Prot No.:	P05067
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	APP Monoclonal Antibody detects endogenous levels of APP protein.
Storage and Stability:	-20°C/1 year