

Antibodies Product Datasheet

Product Name:	ASK 1 Monoclonal Antibody
Catalog Number:	N0049
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCFCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human ASK 1 expressed in E. Coli.
Concentration:	
Gene Name:	MAP3K5
Protein Name:	Mitogen-activated protein kinase kinase kinase 5
Other Name:	MAP3K5; ASK1; MAPKKK5; MEKK5; Mitogen-activated protein kinase kinase kinase 5; Apoptosis signal-regulating kinase 1; ASK-1; MAPK/ERK kinase kinase 5; MEK kinase 5; MEKK 5
Human Swiss Prot No.:	Q99683
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	ASK 1 Monoclonal Antibody detects endogenous levels of ASK 1 protein.
Storage and Stability:	-20°C/1 year