

Antibodies Product Datasheet

Product Name:	ATF-2 Monoclonal Antibody
Catalog Number:	N0051
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF (paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human ATF-2 expressed in E. Coli.
Concentration:	
Gene Name:	ATF2
Protein Name:	Cyclic AMP-dependent transcription factor ATF-2
Other Name:	ATF2; CREB2; CREBP1; Cyclic AMP-dependent transcription factor ATF-2; cAMP-dependent transcription factor ATF-2; Activating transcription factor 2; Cyclic AMP-responsive element-binding protein 2; CREB-2; cAMP-responsive element-binding pro
Human Swiss Prot No.:	P15336
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	ATF-2 Monoclonal Antibody detects endogenous levels of ATF-2 protein.
Storage and Stability:	-20°C/1 year