

Antibodies Product Datasheet

Product Name:	Axl Monoclonal Antibody
Catalog Number:	N0054
Reactivity:	Human
Applications:	WB,IF/ICC,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant extracellular fragment of human Axl fused with hlgGfc tag expressed in HEK293 cell line.
Concentration:	
Gene Name:	AXL
Protein Name:	Tyrosine-protein kinase receptor UFO
Other Name:	AXL; UFO; Tyrosine-protein kinase receptor UFO; AXL oncogene
Human Swiss Prot No.:	P30530
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Axl Monoclonal Antibody detects endogenous levels of Axl protein.
Storage and Stability:	-20°C/1 year