

Antibodies Product Datasheet

Product Name:	B23 Monoclonal Antibody
Catalog Number:	N0055
Reactivity:	Human, Monkey
Applications:	WB, IHC-p, IF/ICELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human B23 expressed in E. Coli.
Concentration:	
Gene Name:	NPM1
Protein Name:	Nucleophosmin
Other Name:	NPM1; NPM; Nucleophosmin; NPM; Nucleolar phosphoprotein B23; Nucleolar protein NO38; Numatrin
Human Swiss Prot No.:	P06748
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	B23 Monoclonal Antibody detects endogenous levels of B23 protein.
Storage and Stability:	-20°C/1 year