

Antibodies Product Datasheet

Product Name:	Bmi-1 Monoclonal Antibody
Catalog Number:	N0071
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCFCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Bmi-1 expressed in E. Coli.
Concentration:	
Gene Name:	BMI1
Protein Name:	Polycomb complex protein BMI-1
Other Name:	BMI1; PCGF4; RNF51; Polycomb complex protein BMI-1; Polycomb group RING finger protein 4; RING finger protein 51
Human Swiss Prot No.:	P35226
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Bmi-1 Monoclonal Antibody detects endogenous levels of Bmi-1 protein.
Storage and Stability:	-20°C/1 year