

Antibodies Product Datasheet

Product Name:	BMPR-II Monoclonal Antibody
Catalog Number:	N0074
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IHC-p, IF/ICELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human BMPR-II expressed in E. Coli.
Concentration:	
Gene Name:	BMPR2
Protein Name:	Bone morphogenetic protein receptor type-2
Other Name:	BMPR2; PPH1; Bone morphogenetic protein receptor type-2; BMP type-2 receptor; BMPR-2; Bone morphogenetic protein receptor type II; BMP type II receptor; BMPR-II
Human Swiss Prot No.:	Q13873
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	BMPR-II Monoclonal Antibody detects endogenous levels of BMPR-II protein.
Storage and Stability:	-20°C/1 year