

Antibodies Product Datasheet

Product Name:	C/EBP β Monoclonal Antibody
Catalog Number:	N0083
Reactivity:	Human
Applications:	WB,FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human C/EBP β expressed in E. Coli.
Concentration:	
Gene Name:	CEBPB
Protein Name:	CCAAT/enhancer-binding protein beta
Other Name:	CEBPB; LAP; TCF5; PP9092; CCAAT/enhancer-binding protein beta; C/EBP beta; Liver activator protein; Nuclear factor NF-IL6; Transcription factor 5; TCF-5
Human Swiss Prot No.:	P17676
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	C/EBP β Monoclonal Antibody detects endogenous levels of C/EBP β protein.
Storage and Stability:	-20°C/1 year