

Antibodies Product Datasheet

Product Name:	CD166 Monoclonal Antibody
Catalog Number:	N0104
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of CD166 (aa405-524) expressed in E. Coli.
Concentration:	
Gene Name:	ALCAM
Protein Name:	CD166 antigen
Other Name:	ALCAM; MEMD; CD166 antigen; Activated leukocyte cell adhesion molecule; CD antigen CD166
Human Swiss Prot No.:	Q13740
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	CD166 Monoclonal Antibody detects endogenous levels of CD166 protein.
Storage and Stability:	-20°C/1 year