

Antibodies Product Datasheet

Product Name:	CD19 Monoclonal Antibody
Catalog Number:	N0107
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human CD19 expressed in E. Coli.
Concentration:	1 mg/ml
Gene Name:	CD19
Protein Name:	B-lymphocyte antigen CD19
Other Name:	CD19; B-lymphocyte antigen CD19; B-lymphocyte surface antigen B4; Differentiation antigen CD19; T-cell surface antigen Leu-12; CD antigen CD19
Human Swiss Prot No.:	P15391
Formulation:	Purified antibody in PBS containing 0.03% sodium azide.
Specificity:	CD19 Monoclonal Antibody detects endogenous levels of CD19 protein.
Storage and Stability:	-20°C/1 year