

Antibodies Product Datasheet

Product Name:	CD276 Monoclonal Antibody
Catalog Number:	N0114
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human CD276 expressed in E. Coli.
Concentration:	
Gene Name:	CD276
Protein Name:	CD276 antigen
Other Name:	CD276; B7H3; PSEC0249; CD276 antigen; 4Ig-B7-H3; B7 homolog 3; B7-H3; Costimulatory molecule; CD antigen CD276
Human Swiss Prot No.:	Q5ZPR3
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	CD276 Monoclonal Antibody detects endogenous levels of CD276 protein.
Storage and Stability:	-20°C/1 year