

Antibodies Product Datasheet

Product Name:	CD31 Monoclonal Antibody
Catalog Number:	N0117
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human CD31 expressed in E. Coli.
Concentration:	
Gene Name:	PDGFRA
Protein Name:	Alpha-type platelet-derived growth factor receptor
Other Name:	PECAM1; Platelet endothelial cell adhesion molecule; PECAM-1; EndoCAM; GPIIA'; PECA1; CD antigen CD31
Human Swiss Prot No.:	P16284
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	CD31 Monoclonal Antibody detects endogenous levels of CD31 protein.
Storage and Stability:	-20°C/1 year