

Antibodies Product Datasheet

Product Name:	CD4 Monoclonal Antibody
Catalog Number:	N0123
Reactivity:	Human
Applications:	FCM,ELISA
Dilutions:	Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human CD4 expressed in E. Coli.
Concentration:	1 mg/ml
Gene Name:	CD4
Protein Name:	T-cell surface glycoprotein CD4
Other Name:	CD4; T-cell surface glycoprotein CD4; T-cell surface antigen T4/Leu-3; CD antigen CD4
Human Swiss Prot No.:	P01730
Formulation:	Purified antibody in PBS containing 0.03% sodium azide.
Specificity:	CD4 Monoclonal Antibody detects endogenous levels of CD4 protein.
Storage and Stability:	-20°C/1 year