

Antibodies Product Datasheet

Product Name:	CD40 Monoclonal Antibody
Catalog Number:	N0125
Reactivity:	Human
Applications:	IF/ICC,ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of CD40 expressed in E. Coli.
Concentration:	1 mg/ml
Gene Name:	CD40
Protein Name:	Tumor necrosis factor receptor superfamily member 5
Other Name:	CD40; TNFRSF5; Tumor necrosis factor receptor superfamily member 5; B-cell surface antigen CD40; Bp50; CD40L receptor; CDw40; CD antigen CD40
Human Swiss Prot No.:	P25942
Formulation:	Purified recombinant fragment of CD40 expressed in E. Coli.
Specificity:	CD40 Monoclonal Antibody detects endogenous levels of CD40 protein.
Storage and Stability:	-20°C/1 year