

Antibodies Product Datasheet

Product Name:	CD66e Monoclonal Antibody
Catalog Number:	N0129
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human CD66e expressed in E. Coli.
Concentration:	
Gene Name:	CEACAM5
Protein Name:	Carcinoembryonic antigen-related cell adhesion molecule 5
Other Name:	CEACAM5; CEA; Carcinoembryonic antigen-related cell adhesion molecule 5; Carcinoembryonic antigen; CEA; Meconium antigen 100; CD antigen CD66e
Human Swiss Prot No.:	P06731
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	CD66e Monoclonal Antibody detects endogenous levels of CD66e protein.
Storage and Stability:	-20°C/1 year