

Antibodies Product Datasheet

Product Name:	CD69 Monoclonal Antibody
Catalog Number:	N0130
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human CD69 expressed in E. Coli.
Concentration:	
Gene Name:	CD69
Protein Name:	Early activation antigen CD69
Other Name:	CD69; CLEC2C; Early activation antigen CD69; Activation inducer molecule; AIM; BL-AC/P26; C-type lectin domain family 2 member C; EA1; Early T-cell activation antigen p60; GP32/28; Leukocyte surface antigen Leu-23; MLR-3; CD antigen CD69
Human Swiss Prot No.:	Q07108
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	CD69 Monoclonal Antibody detects endogenous levels of CD69 protein.
Storage and Stability:	-20°C/1 year