

Antibodies Product Datasheet

Product Name:	CD80 Monoclonal Antibody
Catalog Number:	N0135
Reactivity:	Human
Applications:	IHC-p,IF/ICC,ELISA
Dilutions:	Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of CD80 expressed in E. Coli.
Concentration:	
Gene Name:	CD80
Protein Name:	T-lymphocyte activation antigen CD80
Other Name:	CD80; CD28LG; CD28LG1; LAB7; T-lymphocyte activation antigen CD80; Activation B7-1 antigen; BB1; CTLA-4 counter-receptor B7.1; B7; CD antigen CD80
Human Swiss Prot No.:	P33681
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	CD80 Monoclonal Antibody detects endogenous levels of CD80 protein.
Storage and Stability:	-20°C/1 year