

Antibodies Product Datasheet

Product Name:	CD94 Monoclonal Antibody
Catalog Number:	N0137
Reactivity:	Human
Applications:	WB,FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human CD94 expressed in E. Coli.
Concentration:	
Gene Name:	KLRD1
Protein Name:	Natural killer cells antigen CD94
Other Name:	KLRD1; CD94; Natural killer cells antigen CD94; KP43; Killer cell lectin-like receptor subfamily D member 1; NK cell receptor; CD antigen CD94
Human Swiss Prot No.:	Q13241
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	CD94 Monoclonal Antibody detects endogenous levels of CD94 protein.
Storage and Stability:	-20°C/1 year