

Antibodies Product Datasheet

Product Name:	CDC27 Monoclonal Antibody
Catalog Number:	N0142
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human CDC27 expressed in E. Coli.
Concentration:	
Gene Name:	CDC27
Protein Name:	Cell division cycle protein 27 homolog
Other Name:	CDC27; ANAPC3; D0S1430E; D17S978E; Cell division cycle protein 27 homolog; Anaphase-promoting complex subunit 3; APC3; CDC27 homolog; CDC27Hs; H-NUC
Human Swiss Prot No.:	P30260
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	CDC27 Monoclonal Antibody detects endogenous levels of CDC27 protein.
Storage and Stability:	-20°C/1 year