

Antibodies Product Datasheet

Product Name:	Chitotriosidase Monoclonal Antibody
Catalog Number:	N0151
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of Chitotriosidase (aa22-137) expressed in E. Coli.
Concentration:	
Gene Name:	CHIT1
Protein Name:	Chitotriosidase-1
Other Name:	CHIT1; Chitotriosidase-1; Chitinase-1
Human Swiss Prot No.:	Q13231
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Chitotriosidase Monoclonal Antibody detects endogenous levels of Chitotriosidase protein.
Storage and Stability:	-20°C/1 year