

Antibodies Product Datasheet

Product Name:	Csk Monoclonal Antibody
Catalog Number:	N0169
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IF/ICC, FCM, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Csk expressed in E. Coli.
Concentration:	
Gene Name:	CSK
Protein Name:	Tyrosine-protein kinase CSK
Other Name:	CSK; Tyrosine-protein kinase CSK; C-Src kinase; Protein-tyrosine kinase CYL
Human Swiss Prot No.:	P41240
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	Csk Monoclonal Antibody detects endogenous levels of Csk protein.
Storage and Stability:	-20°C/1 year