

Antibodies Product Datasheet

Product Name:	c-Yes Monoclonal Antibody
Catalog Number:	N0177
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of c-Yes (aa10-193) expressed in E. Coli.
Concentration:	
Gene Name:	YES1
Protein Name:	Proto-oncogene tyrosine-protein kinase Yes
Other Name:	YES1; YES; Tyrosine-protein kinase Yes; Proto-oncogene c-Yes; p61-Yes
Human Swiss Prot No.:	P07947
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	c-Yes Monoclonal Antibody detects endogenous levels of c-Yes protein.
Storage and Stability:	-20°C/1 year