

Antibodies Product Datasheet

Product Name:	Cytokeratin (Pan) Monoclonal Antibody
Catalog Number:	N0190
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of Cytokeratin 5 expressed in E. Coli.
Concentration:	
Gene Name:	KRT5
Protein Name:	Keratin type II cytoskeletal 5
Other Name:	KRT5; Keratin; type II cytoskeletal 5; 58 kDa cytokeratin; Cytokeratin-5; CK-5; Keratin-5; K5; Type-II keratin Kb5
Human Swiss Prot No.:	P13647
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Cytokeratin (Pan) Monoclonal Antibody detects endogenous levels of Cytokeratin (Pan) protein.
Storage and Stability:	-20°C/1 year