

Antibodies Product Datasheet

Product Name:	DDR2 Monoclonal Antibody
Catalog Number:	N0194
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human DDR2 expressed in E. Coli.
Concentration:	
Gene Name:	DDR2
Protein Name:	Discoidin domain-containing receptor 2
Other Name:	DDR2; NTRKR3; TKT; TYRO10; Discoidin domain-containing receptor 2; Discoidin domain receptor 2; CD167 antigen-like family member B; Discoidin domain-containing receptor tyrosine kinase 2; Neurotrophic tyrosine kinase; receptor-related 3; Re
Human Swiss Prot No.:	Q16832
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	DDR2 Monoclonal Antibody detects endogenous levels of DDR2 protein.
Storage and Stability:	-20°C/1 year