

Antibodies Product Datasheet

Product Name:	Egr-1 Monoclonal Antibody
Catalog Number:	N0211
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Egr-1 (aa282-433) expressed in E. Coli.
Concentration:	
Gene Name:	EGR1
Protein Name:	Early growth response protein 1
Other Name:	EGR1; KROX24; ZNF225; Early growth response protein 1; EGR-1; AT225; Nerve growth factor-induced protein A; NGFI-A; Transcription factor ETR103; Transcription factor Zif268; Zinc finger protein 225; Zinc finger protein Krox-24
Human Swiss Prot No.:	P18146
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Egr-1 Monoclonal Antibody detects endogenous levels of Egr-1 protein.
Storage and Stability:	-20°C/1 year