

Antibodies Product Datasheet

Product Name:	EHMT1 Monoclonal Antibody
Catalog Number:	N0212
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of EHMT1 expressed in E. Coli.
Concentration:	1 mg/ml
Gene Name:	EHMT1
Protein Name:	Histone-lysine N-methyltransferase, H3 lysine-9 specific 5
Other Name:	EHMT1; EUHMTASE1; GLP; KIAA1876; KMT1D; Histone-lysine N-methyltransferase EHMT1; Euchromatic histone-lysine N-methyltransferase 1; Eu-HMTase1; G9a-like protein 1; GLP; GLP1; Histone H3-K9 methyltransferase 5; H3-K9-HMTase 5; Lysine N-methy
Human Swiss Prot No.:	Q9H9B1
Formulation:	Purified antibody in PBS containing 0.03% sodium azide.
Specificity:	EHMT1 Monoclonal Antibody detects endogenous levels of EHMT1 protein.
Storage and Stability:	-20°C/1 year