

Antibodies Product Datasheet

Product Name:	eIF4E Monoclonal Antibody
Catalog Number:	N0213
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCFCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human eIF4E expressed in E. Coli.
Concentration:	
Gene Name:	EIF4E
Protein Name:	Eukaryotic translation initiation factor 4E
Other Name:	EIF4E; EIF4EL1; EIF4F; Eukaryotic translation initiation factor 4E; eIF-4E; eIF4E; eIF-4F 25 kDa subunit; mRNA cap-binding protein
Human Swiss Prot No.:	P06730
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	eIF4E Monoclonal Antibody detects endogenous levels of eIF4E protein.
Storage and Stability:	-20°C/1 year