

Antibodies Product Datasheet

Product Name:	Emt Monoclonal Antibody
Catalog Number:	N0217
Reactivity:	Human
Applications:	WB,IF/ICC,FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Emt expressed in E. Coli.
Concentration:	
Gene Name:	ITK
Protein Name:	Tyrosine-protein kinase ITK/TSK
Other Name:	ITK; EMT; LYK; Tyrosine-protein kinase ITK/TSK; Interleukin-2-inducible T-cell kinase; IL-2-inducible T-cell kinase; Kinase EMT; T-cell-specific kinase; Tyrosine-protein kinase Lyk
Human Swiss Prot No.:	Q08881
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Emt Monoclonal Antibody detects endogenous levels of Emt protein.
Storage and Stability:	-20°C/1 year