

Antibodies Product Datasheet

Product Name:	EP-CAM Monoclonal Antibody
Catalog Number:	N0218
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human EP-CAM expressed in E. Coli.
Concentration:	
Gene Name:	EPCAM
Protein Name:	Epithelial cell adhesion molecule
Other Name:	EPCAM; GA733-2; M1S2; M4S1; MIC18; TACSTD1; TROP1; Epithelial cell adhesion molecule; Ep-CAM; Adenocarcinoma-associated antigen; Cell surface glycoprotein Trop-1; Epithelial cell surface antigen; Epithelial glycoprotein; EGP; Epithelial gly
Human Swiss Prot No.:	P16422
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	EP-CAM Monoclonal Antibody detects endogenous levels of EP-CAM protein.
Storage and Stability:	-20°C/1 year